Platelet Aggregation Testing Hospital Draw Guide New Orleans, LA Area Lab



8721 Oak Street New Orleans, LA 70118

(504) 866-7090

Please call our lab at (504) 866-7090 before drawing the patient.

Samples <u>must be received within 2 hrs of draw</u> to complete all testing.

Specimen requirements vary by test (call for pediatric draws):

TEST	NUMBER OF TUBES	SAMPLE TYPE	DAYS PERFORMED
Platelet Aggregation Study	6	3 mL – 3.2% Sodium Citrated, Whole Blood, Room Temperature	By appointment only: Monday - Friday 9am-2pm
Aspirin Sensitivity (ASA)*	4		
Plavix® Sensitivity*	4		
Ristocetin-induced Platelet Aggregation (RIPA)	4		

If ordering the ASA and Plavix Sensitivity together, only 4 x 3mL tubes are required for testing.

Patient Preparation

Fasting: 6 hours

Additional Requirements:

No aspirin** or fish oil for 14 days prior to draw No NSAIDs (e.g. ibuprofen, naproxen, etc.) or antihistamines for 72 hours prior to draw

*Please do not discontinue any prescribed medications without consulting the ordering

provider. Patients may remain on aspirin or Plavix if Aspirin Sensitivity or Plavix Sensitivity testing is ordered.

Instructions:

1. Complete an order form and if available include a list of current medications. 2. At the scheduled appointment, draw the required number of tubes. Invert samples gently to mix. 3. Label the tubes with patient first and last name, DOB, draw date and draw time. 4. No other processing is required. Send immediately at room temperature.

Machaon Diagnostics is a CAP and CLIA-accredited laboratory offering specialized clinical testing and rare disease contract research organization (CRO) services.